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Gants d'examen non stériles – Usage unique

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EN – Instructions for use

Non-sterile examination gloves – Single use



Use: These gloves are mainly designed to protect the hands against microorganisms (bacteria, moulds, viruses) and chemical hazards. Gloves are used on healthy skin. These gloves are for professional use only.

Indication: Protection of the user and the patient when performing care.

Composition: Some gloves may contain components that represent a potential cause of known allergies in some sensitive individuals who may develop contact irritation and/or an allergy. In case of an allergic reaction, it is necessary to consult a doctor. For more information, do not hesitate to contact Euromedis.

Technical characteristics:

- Class I medical device, Category III PPE
- Single-use, non-sterile device.

 PPE 2777	<ul style="list-style-type: none"> - Class I medical device (MD), in accordance with European Regulation (EU) 2017/745. - Category III personal protective equipment (PPE), in accordance with European Regulation (EU) 2016/425. The conformity assessment was carried out by the notified body SATRA Europe (No 2777) 																
	EN ISO 21420:2020 Protective gloves - General requirements and test methods. Read the instructions before use.																
EN ISO 374-1/type A  UVWXYZ <hr style="border-top: 1px dashed black;"/> EN ISO 374-1/type B  XYZ <hr style="border-top: 1px dashed black;"/> EN ISO 374-1/type C 	EN ISO 374-1: 2016+A1: 2018 / Type A, B or C: Protective gloves against dangerous chemicals and microorganisms Type A: more than 30 minutes for at least six test chemicals present in the standard. “UVWXYZ” represent the codes for test chemicals. Type B: more than 30 minutes for at least three test chemicals present in the standard. “XYZ” represent the codes for test chemicals. Type A: more than 10 minutes for at least one test chemical present in the standard. Refer to the label of each glove to find out its type. Permeation levels are based on the following breakthrough times: <table border="1" style="margin: 10px auto; border-collapse: collapse; text-align: center;"> <thead> <tr> <th style="padding: 5px;">Permeation performance level</th> <th style="padding: 5px;">1</th> <th style="padding: 5px;">2</th> <th style="padding: 5px;">3</th> <th style="padding: 5px;">4</th> <th style="padding: 5px;">5</th> <th style="padding: 5px;">6</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Measured breakthrough time (min)</td> <td style="padding: 5px;">> 10</td> <td style="padding: 5px;">> 30</td> <td style="padding: 5px;">> 60</td> <td style="padding: 5px;">> 120</td> <td style="padding: 5px;">> 240</td> <td style="padding: 5px;">> 480</td> </tr> </tbody> </table>			Permeation performance level	1	2	3	4	5	6	Measured breakthrough time (min)	> 10	> 30	> 60	> 120	> 240	> 480
Permeation performance level	1	2	3	4	5	6											
Measured breakthrough time (min)	> 10	> 30	> 60	> 120	> 240	> 480											
EN ISO 374-5:2016  VIRUS	EN ISO 374-5: 2016 / VIRUS Protection against bacteria, fungi and viruses																
	Suitable for food contact. Certified according to European regulations 1935/2004 and 10/2011. Refer to the label of each glove.		Lot number														
	Medical device			Latex-free. Refer to the label of each glove.													
	Do not use if the box is damaged			It should be stored between 10°C and 40°C.													
	AQL is 1.5			Keep out of direct sunlight. Store in a dry place.													

	Date of manufacture and country of origin		Expiry date
	Manufacturer: LABORATOIRES EUROMEDIS		Non-sterile medical device
	Single use only		Follow the sorting instructions

EU-TYPE ASSESSMENT CONDUCTED BY:

SATRA Technology Europe^{LTD}
 Bracetown Business Park, Clonee, Co.
 Meath D15YN2P – Ireland

SATRA Technology Europe is also responsible for the continuous compliance of its products.

Warning:

Chemical resistance data are obtained from assessments carried out under laboratory conditions on palm samples only (except in cases where the glove is equal to or greater than 400 mm - where the cuff is also tested), and refer exclusively to the tested chemical substances. Results may vary due to chemical mixtures. EN ISO 374-4 : 2019 The degradation test results indicate the alteration of the puncture resistance properties of the gloves after exposure to the test chemical. This information does not reflect the actual duration of use in the workplace and does not differentiate between pure chemicals and mixtures. Penetration resistance was assessed under laboratory conditions and only applies to the sample tested. Make sure the gloves are suitable for the intended use. Workplace conditions (temperature, abrasion and degradation) may differ from the typical test. During use, gloves are likely to provide less resistance to the chemical substance due to a change in their physical properties. Movements, snags, friction and degradation caused by exposure to the chemical can significantly reduce the actual duration of use. If corrosive chemicals are used, degradation becomes the primary factor in the choice of chemical resistant gloves. Permeation data of chemicals according to EN ISO 374-1: 2016+A1: 2018 and degradation data according to EN ISO 374-4 : 2019 are available on the package. The declarations of conformity are available on the website: www.laboratoires-euromedis.fr. Any complaints regarding the glove should be sent to service.qualite@euromedis.fr

Instructions for use:

Choose the appropriate size. The gloves are ambidextrous: they can be worn on both hands. Insert five fingers into the cuff and pull the cuff on the wrist as far as possible. Make sure the glove fits snugly around the fingers and palm. The cuff should be snug around the wrist. Put on the glove(s) all the way. To achieve the desired protection, this glove can be used with other PPE or protective medical devices (coveralls, gowns, footwear, respiratory protection equipment, etc.).

Special precautions for dressing and undressing may depend on the field of application (e.g. double gloving).

To remove the gloves, grasp the outside of the glove near the wrist; remove the glove from the hand by turning it inside out. Hold it in the opposite gloved hand. Slide a non-gloved finger under the wrist of the remaining glove, taking care not to touch the outside of the glove. Peel off the remaining glove from the inside to create a “bag” containing both gloves. Dispose of them properly (see **Disposal**).

Precautions for use:

- Do not open the box with a sharp object.
- Before use, check the integrity of the glove (perforation, tear). If gloves tear or puncture during use, discard them. If in doubt, do not use them and take a new pair.
- This glove does not protect against mechanical hazards or ionising radiation. Do not wear these gloves in a containment enclosure.
- Avoid exposing gloves to an open flame.
- The combined use of this device with other higher-class PPE cannot extend the scope of the device. The user is solely responsible for choosing the right combination and adaptation of the glove with other equipment as well as the period during which the glove can be worn, depending on the task to be performed, the level of protection required and the comfort of use.
- Not all food contact gloves are suitable for all foods. With certain types of foods, some gloves may demonstrate excessive migration. For the applicable restrictions and specific foods for which gloves can be used, please contact Euromedis.

- Any serious incident in connection with the use of the device must be reported to the manufacturer and the competent authority of the EU member state where the user resides.
- Latex gloves are contraindicated for people who are allergic to latex. Refer to the label to check the presence or absence of latex.

Storage and duration of use

Keep the gloves in their original packaging, away from humidity, sunlight and light. Keep gloves within the temperature range stated on the package. Do not use gloves beyond the expiry date stated on the package.

Disposal

Depending on the field of application, the user must comply with the national regulations in force for the device and its packaging.

Gloves are for single use only and are not designed to be washed or reused. They must be disposed of after use.

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

SATRA Technology Europe^{LTD}
 Bracetown business Park, Clonee, Co.
 Meath D15YN2P – Irlanda

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